

IL13 Rabbit Pab - P35225

Item Number: CM0022060



Introduction

Rabbit polyclonal antibody against human IL13, validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Suitable for studying allergic inflammation and Th2 immune responses.

[Learn More](#)

Parameter	Value
Product Name	IL13 Rabbit pAb
Remarks/Alias	P600; IL-13; IL13
Species	Human
GeneID (human)	3596
GeneID	3596
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-146 of human IL13 (NP_002179.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35225
Protein Weight	16kDa
Shipping	Ice bag