


IL17A Rabbit Monoclonal Antibody For Ihc-P And Elisa - Q16552

Item Number: CM0002197



Introduction

High-quality rabbit monoclonal antibody against human IL17A (CTLA-8, 18kDa), validated for IHC-P and ELISA. Recognizes human IL17A. Ideal for studying inflammatory responses, antimicrobial defense, and Th17-mediated immunity.

[Learn More](#)

Parameter	Value
Product Name	IL17A Rabbit mAb
Remarks/Alias	IL17; CTLA8; IL-17; ILA17; CTLA-8; IL-17A; IL17A
Species	Human
Gene ID	3605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-155 of human IL17A (NP_002181.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q16552
Protein Weight	18kDa
Shipping	Ice bag