


## IL17A Rabbit Polyclonal Antibody - Q16552

Item Number: CM0016457



### Introduction

IL17A Rabbit Polyclonal Antibody for WB, IF-P, IHC-P, ELISA applications. Detects human, mouse, rat IL-17A. Ideal for studying Th17-mediated inflammation, antimicrobial defense, and epithelial barrier integrity.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | IL17A Rabbit pAb                                                                                                   |
| Remarks/Alias    | IL17; CTLA8; IL-17; ILA17; CTLA-8; IL-17A; IL17A                                                                   |
| Species          | Human                                                                                                              |
| Gene ID          | 3605                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 24-155 of human IL17A (NP_002181.1). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q16552                                                                                                             |
| Protein Weight   | 18kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |