

IL2 Rabbit Monoclonal Antibody For Wb, Elisa - P60568

Item Number: CM0007395



Introduction

Rabbit monoclonal antibody targeting human Interleukin-2 (IL2) for Western blot and ELISA. Suitable for studying T-cell growth, immune tolerance, and immunotherapy research.

[Learn More](#)

Parameter	Value
Product Name	IL2 Rabbit mAb
Remarks/Alias	IL-2; TCGF; lymphokine; IL2
Species	Human
Gene ID (Human)	3558
Gene ID	3558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-153 of human IL2 (NP_000577.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P60568
Protein Weight	17kDa
Shipping	Ice bag