


IL21R Rabbit Pab - Q9Hbe5

Item Number: CM0026737



Introduction

IL21R Rabbit pAb is a high-quality polyclonal antibody for WB and ELISA, detecting human and mouse interleukin-21 receptor (CD360). Ideal for lymphoid biology and disease research.

[Learn More](#)

Parameter	Value
Product Name	IL21R Rabbit pAb
Remarks/Alias	NILR, CD360, IMD56, IL21R
Species	Human
Gene ID (Human)	50615
Gene ID	50615
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-240 of human IL21R (NP_851565.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HBE5
Protein Weight	59kDa
Shipping	Ice bag