


IL4 Rabbit Polyclonal Antibody For Wb, Elisa - P05112

Item Number: CM0017626



Introduction

High-quality rabbit polyclonal antibody targeting human IL-4 (Interleukin-4), validated for WB and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying immune regulation and inflammation.

[Learn More](#)

Parameter	Value
Product Name	IL4 Rabbit pAb
Remarks/Alias	BSF1; IL-4; BCGF1; BSF-1; BCGF-1; IL4
Species	Human
GeneID	3565 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human IL4 (NP_000580.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05112
Protein Weight	17kDa
Shipping	Ice bag