


IL4 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P05112

Item Number: CM0026411



Introduction

Polyclonal antibody targeting human IL4, validated for WB, IHC-P, IF-P, ELISA. Ideal for cytokine research, B-cell activation, and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	IL4 Rabbit pAb
Remarks/Alias	BSF1; IL-4; BCGF1; BSF-1; BCGF-1; IL4
Species	Human
GeneID (Human)	3565
GeneID	3565
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-153 of human IL4 (NP_000580.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P05112
Protein Weight	17kDa
Shipping	Ice bag