

IL4 Rabbit Polymab® For Wb, Elisa - P05112 / P07750

Item Number: CM0006937



Introduction

IL4 Rabbit PolymAb® is a rabbit monoclonal antibody validated for WB and ELISA. It recognizes human and mouse IL4, with a predicted weight of 17 kDa. Ideal for studies on cytokine signaling, immune regulation, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	IL4 Rabbit PolymAb®
Remarks/Alias	BSF1; IL-4; BCGF1; BSF-1; BCGF-1
Species	Human/Mouse
GeneID (Human)	3565
GeneID	3565/16189
Immunogen	Recombinant protein of human/mouse IL4.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05112/P07750
Protein Weight	17kDa
Shipping	Ice bag