


Inmt Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95050

Item Number: CM0014835



Introduction

Rabbit polyclonal antibody to indolethylamine N-methyltransferase (INMT) for WB, IF/ICC, ELISA. Suitable for human and mouse samples, detecting a 29 kDa band. Ideal for methylation and detoxification studies.

[Learn More](#)

Parameter	Value
Product Name	INMT Rabbit pAb
Remarks/Alias	TEMT; INMT
Species	Human
Gene ID (Human)	11185
Gene ID	11185
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-263 of human INMT (NP_006765.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95050
Protein Weight	29kDa
Shipping	Ice bag