

Inos Rabbit Pab - P35228

Item Number: CM0013339



Introduction

Rabbit polyclonal antibody targeting human iNOS (NOS2) for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat. 131kDa. Ideal for inflammation and NO signaling studies.

[Learn More](#)

Parameter	Value
Product Name	iNOS Rabbit pAb
Remarks/Alias	NOS; INOS; NOS2A; HEP-NOS; iNOS
Species	Human
Gene ID (Human)	4843
Gene ID	4843
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 948-1153 of human iNOS (NP_000616.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35228
Protein Weight	131kDa
Shipping	Ice bag