

Insulin Rabbit Polyclonal Antibody For Wb, Ihc-P, If-P, Elisa - P01308

Item Number: CM0019526



Introduction

High-quality rabbit polyclonal antibody targeting human insulin (SWISS P01308), suitable for Western blot, IHC-P, IF-P, and ELISA applications. Cross-reacts with mouse and rat, supporting metabolic disorder research.

[Learn More](#)

Parameter	Value
Product Name	Insulin Rabbit pAb
Remarks/Alias	IDDM; ILPR; IRDN; IDDM1; IDDM2; PNDM4; MODY10; Insulin
Species	Human
Gene ID (Human)	3630
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human INS (NP_000198.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01308
Protein Weight	12 kDa
Shipping	Ice bag