

Integrin Alpha V (Itgav/Cd51) Rabbit Mab For Wb, Ihc-P, Fc (Intra), Elisa - P06756

Item Number: CM0003302



Introduction

Integrin alpha V (ITGAV/CD51) rabbit monoclonal antibody for human, mouse, rat. Validated for Western blot, IHC-P, flow cytometry, ELISA. Unconjugated, ~116kDa target, synthetic peptide immunogen.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha V (ITGAV/CD51) Rabbit mAb
Remarks/Alias	CD51; MSK8; VNRA; VTNR; Integrin alpha V (ITGAV/CD51)
Species	Human
GeneID (Human)	3685
GeneID	3685
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 949-1048 of human Integrin alpha V (ITGAV/CD51) (NP_002201.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06756
Protein Weight	116kDa
Shipping	Ice bag