

Integrin Linked Ilk Rabbit Monoclonal Antibody For Wb, Elisa - Q13418

Item Number: CM0007837



Introduction

Rabbit monoclonal antibody targeting human Integrin linked ILK (ILK, P59, ILK-1, ILK-2). Validated for WB and ELISA; cross-reacts with human, mouse, rat. Detects a 51 kDa scaffold protein critical for integrin signaling and focal adhesion dynamics.

[Learn More](#)

Parameter	Value
Product Name	Integrin linked ILK Rabbit mAb
Remarks/Alias	P59; ILK-1; ILK-2; p59ILK; HEL-S-28; Integrin linked ILK
Species	Human
Gene ID (Human)	3611
Gene ID	3611
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 350-452 of human Integrin linked ILK (Q13418).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13418
Protein Weight	51kDa
Shipping	Ice bag