

Ip3R Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q14643

Item Number: CM0012402

Introduction

Rabbit monoclonal antibody against human IP3R (ITPR1). Validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Targets calcium channel involved in intracellular calcium signaling and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	IP3R Rabbit mAb
Remarks/Alias	ACV; CLA4; IP3R; IP3R1; SCA15; SCA16; SCA29; INSP3R1; PPP1R94
Species	Human
GeneID (Human)	3708
GeneID	3708
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 2200-2300 of human IP3R (Q14643).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14643
Protein Weight	314kDa
Shipping	Ice bag