

Ipo11 Rabbit Pab For Wb And Elisa - Q9Ui26

Item Number: CM0015976

Introduction

IPO11 Rabbit Polyclonal Antibody validated for WB and ELISA in human and mouse. Targets Importin-11, a nuclear transport receptor involved in Ran-dependent protein import. For research use.

[Learn More](#)

Parameter	Value
Product Name	IPO11 Rabbit pAb
Remarks/Alias	RanBP11; IPO11
Species	Human
GeneID (Human)	51194
GeneID	51194
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 776-975 of human IPO11 (NP_057422.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Ui26
Protein Weight	113kDa
Shipping	Ice bag