


Irak2 Rabbit Polyclonal Antibody For Wb, Elisa - O43187

Item Number: CM0028839



Introduction

IRAK2 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat IRAK2. Ideal for IL-1 signaling and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	IRAK2 Rabbit pAb
Remarks/Alias	IRAK-2; IRAK2
Species	Human
GeneID (human)	3656
GeneID	3656
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human IRAK2 (NP_001561.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43187
Protein Weight	69kDa
Shipping	Ice bag