


Irak3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y616

Item Number: CM0013864



Introduction

IRAK3 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects endogenous IRAK3 in human, mouse, and rat samples. Ideal for studying IL-1 receptor and Toll-like receptor signaling pathways. - Q9Y616

[Learn More](#)

Parameter	Value
Product Name	IRAK3 Rabbit pAb
Remarks/Alias	ASRT5; IRAKM; IRAK3
Species	Human
GeneID (Human)	11213
GeneID	11213
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 367-596 of human IRAK3 (NP_009130.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y616
Protein Weight	68kDa
Shipping	Ice bag