


Isg15 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05161

Item Number: CM0013107



Introduction

ISG15 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human and rat samples. Ubiquitin-like protein ISG15 plays key roles in innate antiviral immunity and ISGylation. Ideal for IVD and research.

[Learn More](#)

Parameter	Value
Product Name	ISG15 Rabbit pAb
Remarks/Alias	G1P2; IP17; UCRP; IFI15; IMD38; hUCRP; ISG15
Species	Human
Gene ID (Human)	9636
Gene ID	9636
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human ISG15 (NP_005092.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P05161
Protein Weight	18kDa
Shipping	Ice bag