

Itgb3 Rabbit Polyclonal Antibody For Wb, Fc, Elisa - P05106

Item Number: CM0024490



Introduction

ITGB3 Rabbit Polyclonal Antibody validated for WB, FC, ELISA. Reacts with human, mouse, rat. Integrin beta-3 (CD61) is essential for platelet aggregation and cell adhesion. Suitable for diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	ITGB3 Rabbit pAb
Remarks/Alias	ITGB3; BDPLT16; BDPLT2; CD61; GP3A; GPIIIa; GT; integrin beta-3
Species	Human
GeneID (Human)	3690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-718 of human ITGB3 (NP_000203.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05106
Protein Weight	86kDa/87kDa
Shipping	Ice bag