

Itpk1 Rabbit Pab - Q13572

Item Number: CM0025819



Introduction

High-quality rabbit polyclonal antibody against human ITPK1, validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying inositol phosphate metabolism, apoptosis, and necroptosis.

[Learn More](#)

Parameter	Value
Product Name	ITPK1 Rabbit pAb
Remarks/Alias	ITRPK1; ITPK1
Species	Human
Gene ID (Human)	3705
Gene ID	3705
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 185-414 of human ITPK1 (NP_055031.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13572
Protein Weight	46kDa
Shipping	Ice bag