

Jab1/Csn5/Cops5 Rabbit Polyclonal Antibody For Wb And Elisa - Q92905

Item Number: CM0022270



Introduction

Rabbit polyclonal antibody against human JAB1/CSN5/COPS5 (COP9 signalosome subunit 5). Validated for WB and ELISA; cross-reacts with human, mouse, and rat. Enables study of ubiquitin conjugation and deneddylation pathways.

[Learn More](#)

Parameter	Value
Product Name	JAB1/CSN5/COPS5 Rabbit pAb
Remarks/Alias	CSN5; JAB1; SGN5; MOV-34; JAB1/CSN5/COPS5
Species	Human
Gene ID (Human)	10987
Gene ID	10987
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human JAB1/CSN5/COPS5 (NP_006828.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92905
Protein Weight	38kDa
Shipping	Ice bag