

Jak2 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - O60674

Item Number: CM0016063



Introduction

High-quality rabbit polyclonal antibody against human JAK2, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for JAK/STAT signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	JAK2 Rabbit pAb
Remarks/Alias	JTK10; JAK2
Species	Human
GeneID (Human)	3717
GeneID	3717
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 933-1132 of human JAK2 (NP_004963.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60674
Protein Weight	131kDa
Shipping	Ice bag