


## Jnk3 Rabbit Monoclonal Antibody For Wb And Elisa - P53779

Item Number: CM0003495



### Introduction

Rabbit monoclonal antibody targeting JNK3/MAPK10, validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for studies on stress-activated signaling, neuronal apoptosis, and circadian rhythms.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | JNK3 Rabbit mAb                                                                                      |
| Remarks/Alias    | JNK3; JNK3A; PRKM10; SAPK1b; p493F12; p54bSAPK                                                       |
| Species          | Human                                                                                                |
| Gene ID (Human)  | 5602                                                                                                 |
| Gene ID          | 5602                                                                                                 |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 350-450 of human MAPK10 (P53779). |
| Source           | Rabbit                                                                                               |
| Category         | Monoclonal Antibodies                                                                                |
| Application      | WB, ELISA                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                    |
| SWISS            | P53779                                                                                               |
| Protein Weight   | 44kDa                                                                                                |
| Shipping         | Ice bag                                                                                              |