

K48-Linkage Specific Polyubiquitin Rabbit Polyclonal Antibody For Wb, Dot Blot, Elisa

Item Number: CM0000033

Introduction

K48-linkage specific polyubiquitin rabbit polyclonal antibody suitable for Western blot, dot blot, and ELISA. Predicted to react with human, mouse, rat, and a wide range of species. Useful for ubiquitin-proteasome pathway studies.

[Learn More](#)

Parameter	Value
Product Name	K48-linkage Specific Polyubiquitin Rabbit pAb
Gene ID (Human)	7314
Immunogen	Synthetic peptide corresponding to a sequence containing polyubiquitinated protein (K48 linkage).
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag