

K48-Linkage Specific Ubiquitin Rabbit Mab - P0Cg47

Item Number: CM0005308



Introduction

This rabbit monoclonal antibody specifically detects K48-linked polyubiquitin chains in human samples. Suitable for WB, IF/ICC, and ELISA applications. Target protein weight ~26kDa (UniProt P0CG47).

[Learn More](#)

Parameter	Value
Product Name	K48-linkage Specific Ubiquitin Rabbit mAb
Remarks/Alias	HEL-S-50; K48-linkage Specific Ubiquitin
Species	Human
Gene ID	7314
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence containing Ubiquitin protein (K48 linkage)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P0CG47
Protein Weight	26kDa
Shipping	Ice bag