

K63-Linkage Specific Polyubiquitin Rabbit Polyclonal Antibody For Dot Blot And Elisa Applications

Item Number: CM0000059



Introduction

This rabbit polyclonal antibody specifically targets K63-linked polyubiquitinated proteins, suitable for dot blot and ELISA. It is predicted to cross-react widely with human, mouse, rat, and other species.

[Learn More](#)

Parameter	Value
Product Name	K63-linkage Specific Polyubiquitin Rabbit pAb
GeneID (Human)	7314
Immunogen	A synthetic peptide corresponding to a sequence containing polyubiquitinated protein (K63 linkage)
Source	Rabbit
Category	Other Modified Antibodies
Application	DB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag