

Kbtbd3 Rabbit Pab - Np_081238.2

Item Number: CM0023149

Introduction

KBTBD3 Rabbit pAb is a polyclonal antibody raised in rabbit against a recombinant fragment of mouse KBTBD3. Validated for WB and ELISA, with cross-reactivity to human. Suitable for ubiquitination and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	KBTBD3 Rabbit pAb
Remarks/Alias	Bklhd3; Z200003A07Rik; KBTBD3
Species	Mouse
Gene ID (Human)	143879
Gene ID (Mouse)	69149
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 410-607 of mouse KBTBD3 (NP_081238.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
Protein Weight	69kDa
Shipping	Ice bag