

Kbtbd7 Rabbit Pab - Q8Wvz9

Item Number: CM0027819

Introduction

KBTBD7 Rabbit pAb is a polyclonal antibody for WB and ELISA, targeting human KBTBD7 with cross-reactivity to mouse and rat. Investigates RAC1 signaling, ubiquitination, and cytoskeletal regulation.

[Learn More](#)

Parameter	Value
Product Name	KBTBD7 Rabbit pAb
Remarks/Alias	KBTBD7
Species	Human
Gene ID (Human)	84078
Gene ID	84078
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 475-684 of human KBTBD7 (NP_115514.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WVZ9
Protein Weight	77kDa
Shipping	Ice bag