

Kcnj3 Rabbit Pab For Wb, Elisa - P48549

Item Number: CM0027722

Introduction

Rabbit polyclonal antibody targeting human KCNJ3 (GIRK1), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying G protein-activated inward rectifier potassium channels.

[Learn More](#)

Parameter	Value
Product Name	KCNJ3 Rabbit pAb
Remarks/Alias	KGA; GIRK1; KIR3.1; KCNJ3
Species	Human
GeneID (human)	3760
GeneID	3760
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human KCNJ3 (NP_002230.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48549
Protein Weight	57kDa
Shipping	Ice bag