


Kcnk9 Polyclonal Antibody (Rabbit) For Wb, Elisa - Q9Npc2

Item Number: CM0029373



Introduction

High-quality KCNK9 rabbit polyclonal antibody for Western blot and ELISA. Reacts with human, mouse, and rat KCNK9, a pH-sensitive potassium channel involved in neuronal excitability and hormone regulation. Ideal for neuroscience and ion channel studies.

[Learn More](#)

Parameter	Value
Product Name	KCNK9 Rabbit pAb
Remarks/Alias	KT3.2; TASK3; BIBARS; K2p9.1; TASK-3; TASK32; KCNK9
Species	Human
GeneID (Human)	51305
GeneID	51305
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KCNK9 (NP_001269463.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPC2
Protein Weight	42kDa
Shipping	Ice bag