


Kcnmb2 Rabbit Pab For Wb, Elisa - Q9Y691

Item Number: CM0016967



Introduction

KCNMB2 Rabbit Polyclonal Antibody validated for WB and ELISA, reactive with human, mouse, and rat. Targets BK channel beta-2 subunit, modulating KCNMA1 potassium channel inactivation. Suitable for ion channel research in neurological and renal tissues.

[Learn More](#)

Parameter	Value
Product Name	KCNMB2 Rabbit pAb
Remarks/Alias	KCNMB2
Species	Human
Gene ID (Human)	10242
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 68-194 of human KCNMB2 (NP_005823.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y691
Protein Weight	27kDa
Shipping	Ice bag