


Kcnu1 Rabbit Polyclonal Antibody For Wb And Elisa - A8Myu2

Item Number: CM0014840



Introduction

Rabbit polyclonal antibody against human KCNU1 (Slo3/KCa5). Validated for WB, ELISA; cross-reacts with rat. Targets testis-specific potassium channel critical for sperm capacitation and fertility.

[Learn More](#)

Parameter	Value
Product Name	KCNU1 Rabbit pAb
Remarks/Alias	KCa5; Slo3; KCNMC1; KCa5.1; Kcnma3; SPGF79; KCNU1
Species	Human
GeneID (Human)	157855
GeneID	157855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 920-1149 of human KCNU1 (NP_001027006.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	A8MYU2
Protein Weight	130kDa
Shipping	Ice bag