

[Kd Validated] Anti-Acsl4 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - O60488

Item Number: CM0007060

Introduction

KD validated recombinant rabbit monoclonal antibody specific for ACSL4. For Western blot, IF/ICC, IHC-P, ELISA. Detects human, mouse, and rat ACSL4. Key tool for ferroptosis and fatty acid metabolism research.

[Learn More](#)

Parameter	Description
Product Name	[KD Validated] ACSL4 Rabbit mAb
Remarks/Alias	ACS4; FAACL4; LACS4; MRX63; MRX68; XLID63
Species	Human
GeneID (Human)	2182
GeneID	2182
Immunogen	Synthetic peptide derived from human ACSL4
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60488
Protein Weight	79 kDa
Shipping	Ice bag