

[Kd Validated] Anti-Acss2 Rabbit Polymab® For Wb, Elisa - Q9Nr19

Item Number: CM0002540

Introduction

Validated anti-ACSS2 rabbit monoclonal antibody for Western blot and ELISA. Recognizes human, mouse, rat ACSS2. Suitable for studies on cytoplasmic acetyl-CoA metabolism, autophagy, and tumorigenesis.

[Learn More](#)

Parameter	Details
Product Name	[KD Validated] ACSS2 Rabbit PolymAb®
Remarks/Alias	ACS, ACSA, ACAS2, ACECS, AceCS1, dj1161H23.1
Species	Human
GeneID (Human)	55902
GeneID	55902
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 50-150 of human ACSS2 (Q9NR19)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR19
Protein Weight	79 kDa/80 kDa
Shipping	Ice bag