


[Kd Validated] Anti-Adar1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P55265

Item Number: CM0011550



Introduction

This [KD Validated] recombinant rabbit monoclonal antibody targets human ADAR1 (P55265, 136 kDa). Validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for RNA editing and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] ADAR1 Rabbit mAb
Remarks/Alias	DSH; AGS6; G1P1; IFI4; P136; ADAR1; DRADA; DSRAD; IFI-4; K88DSRBP
Species	Human
Gene ID (Human)	103
Gene ID	103
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 610-810 of human ADAR1 (NP_001102.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P55265
Protein Weight	136kDa
Shipping	Ice bag