

Kd Validated Anti-Atf1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P18846

Item Number: CM0011458

Introduction

KD validated rabbit monoclonal antibody against human ATF1 (TREB36, EWS-ATF1). Reacts with human, mouse, rat. Suitable for WB, IF/ICC, IHC-P, ELISA. Ideal for studying transcription regulation, cell proliferation, and stress responses.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] ATF1 Rabbit mAb
Remarks/Alias	TREB36; EWS-ATF1; FUS/ATF-1
Species	Human
GeneID (Human)	466
GeneID	466
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18846
Protein Weight	29kDa
Shipping	Ice bag