

[Kd Validated] Anti-Atg7 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - O95352

Item Number: CM0003628

Introduction

KD validated anti-ATG7 rabbit monoclonal antibody for Western blot, IHC-P, IP, and ELISA. Detects human, mouse, and rat ATG7 (~78 kDa). Supports autophagy and mitophagy research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] ATG7 Rabbit mAb
Remarks/Alias	GSA7; [KD Validated] ATG7
Species	Human
Gene ID (Human)	10533
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Atg7(Apg7) (O95352).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95352
Protein Weight	78 kDa
Shipping	Ice bag