

[Kd Validated] Anti-Aurora B Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q96Gd4

Item Number: CM0004060



Introduction

Reliable [KD Validated] anti-Aurora B rabbit monoclonal antibody suitable for WB, IF/ICC, and ELISA applications. Detects Aurora kinase B in human, mouse, and rat samples. Ideal for studies on mitosis, chromosome segregation, and cytokinesis.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Aurora B Rabbit mAb
Remarks/Alias	AIK2; AIM1; ARK2; AurB; IPL1; STK5; AIM-1; ARK-2; STK-1; STK12; PPP1R48; aurkb-sv1; aurkb-sv2; [KD Validated] Aurora B
Species	Human
GeneID (human)	9212
GeneID	9212
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Aurora B (NP_004208.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96GD4
Protein Weight	39kDa
Shipping	Ice bag