

[Kd Validated] Anti-Bcl-2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P10415

Item Number: CM0002068



Introduction

High-specificity [KD Validated] Bcl-2 rabbit monoclonal antibody validated for WB, IHC-P, IP, ELISA. Recognizes human and mouse Bcl-2 with a predicted weight of 26 kDa. Ideal for apoptosis and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Bcl-2 Rabbit mAb
Remarks/Alias	Bcl-2; PPP1R50
Species	Human
GeneID (Human)	596
GeneID	596
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human Bcl-2 (P10415).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10415
Protein Weight	26 kDa
Shipping	Ice bag