


## Kd Validated Anti-Bub1 Monoclonal Antibody For Wb, Elisa - O43683

Item Number: CM0006775



### Introduction

KD-validated rabbit monoclonal antibody against human BUB1, suitable for Western blot and ELISA. Essential for spindle-assembly checkpoint and chromosome alignment during mitosis. Cross-reacts with human. Protein weight ~115-122 KDa.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] BUB1 Rabbit mAb                                                                                   |
| Remarks/Alias    | BUB1A, BUB1L, hBUB1, MCPH30, KD Validated BUB1                                                                   |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 699                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human BUB1. (NP_004327.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | O43683                                                                                                           |
| Protein Weight   | 115 KDa / 120 KDa / 122 KDa                                                                                      |
| Shipping         | Ice bag                                                                                                          |