

[Kd Validated] Anti-Casein Kinase 2 Alpha (Csnk2A1) Rabbit Monoclonal Antibody For Wb, Elisa - P68400

Item Number: CM0003795



Introduction

Anti-Casein Kinase 2 alpha rabbit monoclonal antibody validated for WB and ELISA. Targets human CSNK2A1, a key kinase in cell cycle, apoptosis, transcription, and DNA repair pathways. Research use, not for diagnostics.

[Learn More](#)

Parameter	Specification
Product Name	[KD Validated] Casein Kinase 2 alpha (CSNK2A1) Rabbit mAb
Remarks/Alias	CKII, Cka1, Cka2, CK2A1, OCNDS, [KD Validated] Casein Kinase 2 alpha (CSNK2A1)
Species	Human
GeneID (human)	1457
GeneID	1457
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P68400
Protein Weight	45 kDa
Shipping	Ice bag