

[Kd Validated] Anti-Cat1 Rabbit Polyclonal Antibody For Wb - P30825

Item Number: CM0028158

Introduction

High-affinity cationic amino acid transporter 1 (CAT1/SLC7A1) rabbit polyclonal antibody validated for Western blot and ELISA applications; detects human and mouse CAT1. Suitable for transport and viral receptor studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CAT1 Rabbit pAb
Remarks/Alias	ERR; ATRC1; CAT-1; HCAT1; REC1L
Species	Human
GeneID (Human)	6541
GeneID	6541
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human SLC7A1 (NP_003036.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30825
Protein Weight	68kDa
Shipping	Ice bag