


[Kd Validated] Anti-Ccnt1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - O60563

Item Number: CM0007726



Introduction

Validated rabbit monoclonal antibody against human Cyclin T1 (CCNT1), essential for transcriptional elongation. Suitable for WB, IP, IF/ICC, and ELISA. Ideal for transcription and HIV research applications.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CCNT1 Rabbit mAb
Remarks / Alias	CCNT; CYCT1; HIVE1
Species	Human
Gene ID (Human)	904
Gene ID	904
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 621-720 of human Cyclin T1 (NP_001231.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O60563
Protein Weight	81kDa
Shipping	Ice bag