


[Kd Validated] Anti-Cd98/Slc3A2 Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P08195

Item Number: CM0005096



Introduction

KD Validated rabbit monoclonal antibody against human CD98/SLC3A2 (4F2hc), validated for WB, IF/ICC, FC, ELISA. Essential for amino acid transport, integrin signaling, and viral entry. Ideal for cancer and infection research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CD98/SLC3A2 Rabbit mAb
Remarks/Alias	4F2; CD98; MDU1; 4F2HC; 4T2HC; NACAE; CD98HC; CD98/SLC3A2
Species	Human
GeneID (Human)	6520
GeneID	6520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-630 of human CD98/SLC3A2 (NP_002385.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P08195
Protein Weight	57kDa/61kDa/67kDa/68kDa/71kDa
Shipping	Ice bag