

[Kd Validated] Anti-Cox Iv Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P13073

Item Number: CM0019735



Introduction

This rabbit polyclonal antibody targets human cytochrome c oxidase subunit IV (COX IV), validated for WB, IHC-P, IF/ICC, and ELISA. It cross-reacts with mouse and rat, ideal for studying mitochondrial electron transport and oxidative phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] COX IV Rabbit pAb
Remarks/Alias	COX4; COXIV; COX4-1; COXIV-1; MC4DN16; COX IV-1; [KD Validated] COX IV
Species	Human
GeneID (human)	1327
GeneID	1327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-98 of human COX IV (NP_001852.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13073
Protein Weight	20kDa
Shipping	Ice bag