


[Kd Validated] Anti-Crebbp Rabbit Monoclonal Antibody For Wb/If/Ihc - Q92793

Item Number: CM0009257



Introduction

Reliable rabbit monoclonal CREBBP antibody validated for knockdown applications. Detects human, mouse, and rat CREBBP in WB, IF/ICC, IHC-P, and ELISA. UniProt: Q92793.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CREBBP Rabbit mAb
Remarks/Alias	CBP; RSTS; KAT3A; MKHK1; RSTS1
Species	Human
GeneID (Human)	1387
GeneID	1387
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 2343-2442 of human CREBBP (NP_004371.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92793
Protein Weight	265kDa
Shipping	Ice bag