


[Kd Validated] Anti-Cyb5R3 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00387

Item Number: CM0005014



Introduction

KD-validated rabbit monoclonal antibody targeting human CYB5R3 (DIA1/B5R) for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. 34kDa protein. Ideal for metabolism and erythroid research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CYB5R3 Rabbit mAb
Remarks/Alias	B5R; DIA1
Species	Human
GeneID (Human)	1727
Immunogen	Recombinant protein corresponding to a sequence within amino acids 27-301 of human CYB5R3 (NP_000389.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00387
Protein Weight	34kDa
Shipping	Ice bag