


Kd Validated Anti-Cyclophilin F Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P30405

Item Number: CM0009828



Introduction

KD Validated rabbit monoclonal antibody targeting human Cyclophilin F (PPIF, CypD). Applications: WB, IHC-P, ELISA. Ideal for mitochondrial apoptosis, necrosis, and permeability transition pore research. Molecular weight 22kDa, cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Cyclophilin F Rabbit mAb
Remarks/Alias	CYP3; CypD; CyP-M; Cyp-D
Species	Human
GeneID (Human)	10105
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 108-207 of human Cyclophilin F (NP_005720.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P30405
Protein Weight	22kDa
Shipping	Ice bag