

[Kd Validated] Anti-Dec1 Rabbit Monoclonal Antibody For Wb, Elisa - O14503

Item Number: CM0007310

Introduction

High-quality rabbit monoclonal antibody against human DEC1 (BHLHE40), validated for Western blot and ELISA. Designed for circadian rhythm, transcription, and chondrocyte differentiation research. SWISS: O14503.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] DEC1 Rabbit mAb
Remarks/Alias	DEC1; HLHB2; BHLHB2; Clast5; SHARP2; STRA13; Stra14; SHARP-2
Species	Human
GeneID (human)	8553
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 121-300 of human DEC1 (NP_003661.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14503
Protein Weight	46 kDa
Shipping	Ice bag