


Kd Validated Anti-Dync1H1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14204

Item Number: CM0009089



Introduction

KD validated rabbit monoclonal antibody targeting human DYNC1H1. Suitable for WB, IHC-P, IF/ICC, ELISA applications. Reacts with human, mouse, and rat. Specific to cytoplasmic dynein heavy chain 1.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] DYNC1H1 Rabbit mAb
Remarks/Alias	p22; DHC1; DNCL; DYHC; HL-3; CMT2O; DHC1a; DNCH1; DNECL; Dnchc1; CDCBM13; SMALED1
Species	Human
GeneID (Human)	1778
GeneID	1778
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human DYNC1H1 (NP_001367.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14204
Protein Weight	532kDa
Shipping	Ice bag